





Made In Switzerland

implantswiss
IMPLANT TECHNOLOGY

SMART
CHOICE


OVERVIEW



Switzerland is at the heart of
science, mechanical design and craftsmanship

Switzerland has been one of the most important centers for craftsmanship that requires:
high precision and expertise throughout history.

The products with impeccable mechanical designs, can be used for a lifetime and add value to your life.



Implant**swiss** Dental System owes its **excellent** structure to the inspiration, experience and value from history.

Implantswiss**** is a strong and reliable brand in the dental industry with its innovative solutions, unique surface design and high quality.

Implant**swiss** is devoted to Research and Development through cooperations with universities and crucial investments in clinical research.

Implant**swiss** follows the latest trends and initiates extensive and contemporary projects that may lead the future of implant dentistry.

Implant**swiss** products are manufactured with the state-of-the-arts technologies and high precision.

Implant**swiss** never compromises on its quality and reliability; the success and high quality of Implant**swiss** lies in this perfectionist approach in each step of the manufacturing process.

Each Implant**swiss** product can be easily traced from production to delivery.

Each product package has a barcode label with a lot number created by the ERP software system. All steps from production to stock, delivery and after sales can be monitored with the barcode system; therefore, the follow-up process can be conducted flawlessly.



ø 2.9

ø 3.3

ø 3.7

ø 4.3

ø 4.8

ø 5.5

PLATFORM COLOR CODING

CERTIFICATES



CE

ISO

FDA



BONE LEVEL




Conical connection with hexagonal index

Creates an excellent implant-abutment connection.



Platform Switching

Protects the alveolus




Self Cutting Edge

Unique thread design that does not transfer stress to the bone.



Hygienic Neck Design

Easy to maintain the hygiene



STANDARD LINE

Bone Level implants are designed to provide excellent outcome in any bone quality and in any location.


Implantswiss Bone Level implants allow to design an ideal smile in the aesthetic zone.

	ø2.9	ø3.3	ø3.7	ø4.3	ø4.8	ø5.5
Length	ø2.9	ø3.3	ø3.7	ø4.3	ø4.8	ø5.5
8.0 mm	-	S-BMFSR3308	S-BFHR3708	S-BFHR4308	S-BFHR4808	S-BWFRH5508
10.0 mm	S-BNFRH2910	S-BMFSR3310	S-BFHR3710	S-BFHR4310	S-BFHR4810	S-BWFRH5510
12.0 mm	S-BNFRH2912	S-BMFSR3312	S-BFHR3712	S-BFHR4312	S-BFHR4812	S-BWFRH5512
14.0 mm	S-BNFRH2914	S-BMFSR3314	S-BFHR3714	S-BFHR4314	S-BFHR4814	-


SRA* Surface

Following sandblasting, the surface of the implants is etched using Double Acid Etching technique. Implantswiss S.R.A. surface has a homogenous surface roughness that perfectly mimics the bone morphology.


* Sandblasted Rough Acid-etched Surface



3D-Area Parameters



Detail of the thread bank and representation of the surface shape after removal of a 4th order polynomial



Detail of the thread valley and representation of the surface shape after removal of a 4th order polynomial

Hybrid Design

It is a unique implant design with cylindrical (straight) design on the neck area and conical (tapered) design through the apex of the implant.




Straight form

Tapered form + Straight form

Tapered form

SHORT IMPLANT

In cases of insufficient bone volume, either above the mandibular canal or below the maxillary sinuses, short implants may be an alternative treatment to advanced surgical interventions aiming to regenerate hard tissues.

	ø4.3	ø4.8	ø5.5
			
Length	ø4.3	ø4.8	ø5.5
5.0 mm	S-BMFSR4305	S-BMFSR4805	S-BMFSR5505
6.0 mm	S-BMFSR4306	S-BMFSR4806	S-BMFSR5506

ZYGOMA IMPLANTS

ZYGOMA Implant

	ø4.0
35 mm	S-ZYG35
37.5 mm	S-ZYG37.5
40 mm	S-ZYG40
42.5 mm	S-ZYG42.5
45 mm	S-ZYG45
47.5 mm	S-ZYG47.5
50 mm	S-ZYG50
52.5 mm	S-ZYG52.5
55 mm	S-ZYG55
57.5 mm	S-ZYG57.5
60 mm	S-ZYG60



MULTI UNIT

MULTI UNIT ABUTMENT

Mini



ø4.8

only for ø3.3 implant

0.5 mm	S-BMMUA4805
1.5 mm	S-BMMUA4815
2.5 mm	S-BMMUA4825
3.5 mm	S-BMMUA4835
4.5 mm	S-BMMUA4845

Straight



ø4.8

S-BMUA4805
S-BMUA4815
S-BMUA4825
S-BMUA4835
S-BMUA4845

MULTI UNIT ABUTMENT Angled



ø4.8 (17°)



ø4.8 (35°)

1.0 mm	S-BMUA48171	S-BMUA48352
2.5 mm	S-BMUA48351	S-BMUA48353
3.5 mm	S-BMUA48172	S-BMUA48353



MULTI UNIT ABUTMENT PLASTIC CYLINDER

S-MUAPICYN



MULTI UNIT ABUTMENT PLASTIC CYLINDER

S-MUAPICY



MULTI UNIT ABUTMENT TITANIUM CYLINDER

S-MUATICYN



MULTI UNIT IMPRESSION COPING PICK UP

S-BMUAICP ø5.0



MULTI UNIT IMPRESSION COPING TRANSFER

S-MUATN



MULTI UNIT LAB ANALOG

S-BMUALA



MULTI UNIT ABUTMENT HEALING CAP

S-MUAHCap ø4.8

DIGITAL DENTISTRY

BONE LEVEL

- Couple System
Scan Body

- Ti-Base Scan Body

- Digital Analog

Scan Body for Multi Unit

- Scan Body - Digital Analog



only for ø3.3 implant

S-PMSCB



S-PSB



S-TIMSB



S-TISB



only for ø3.3 implant

S-BMFDLA



S-BFDLA



S-PMUSB



S-BMDLA

COUPLE SYSTEM

BONE LEVEL COUPLE ABUTMENTS



only for ø3.3 implant

	1 mm	2 mm	3 mm	4 mm	5 mm
ø3.7	S-BMCA37107	S-BMCA37207	S-BMCA37307	S-BMCA37407	S-BMCA37507
ø4.5	S-BCA45107	S-BCA45207	S-BCA45307	S-BCA45407	S-BCA45507
ø5.5	S-BCA55107	S-BCA55207	S-BCA55307	S-BCA55407	S-BCA55507
ø6.5	S-BCA65107	S-BCA65207	S-BCA65307	S-BCA65407	S-BCA65507

BONE LEVEL COUPLE ABUTMENTS NO HEX



only for ø3.3 implant

	1 mm	2 mm	3 mm	4 mm	5 mm
ø3.7	S-BMCAN37107	S-BMCAN37207	S-BMCAN37307	S-BMCAN37407	S-BMCAN37507
ø4.5	S-BCAN45107	S-BCAN45207	S-BCAN45307	S-BCAN45407	S-BCAN45507
ø5.5	S-BCAN55107	S-BCAN55207	S-BCAN55307	S-BCAN55407	S-BCAN55507
ø6.5	S-BCAN65107	S-BCAN65207	S-BCAN65307	S-BCAN65407	S-BCAN65507

BONE LEVEL ANGLED ABUTMENTS



only for ø3.3 implant

	1 mm (15°)	2 mm (15°)	4 mm (15°)
ø3.7	S-BMAA37151	S-BMAA37152	S-BMAA37154
ø4.5	S-BAA45152	S-BAA45252	S-BAA45154
ø5.5	S-BAA55152	S-BAA55252	S-BAA55154
ø6.5	S-BAA65152	S-BAA65252	S-BAA65154

BONE LEVEL ANGLED ABUTMENTS NO HEX



only for ø3.3 implant

	1 mm (15°)	2 mm (15°)	4 mm (15°)
ø3.7	S-BMAAN37151	S-BMAAN37152	S-BMAAN37154
ø4.5	S-BAAN45152	S-BAAN45252	S-BAAN45154
ø5.5	S-BAAN55152	S-BAAN55252	S-BAAN55154
ø6.5	S-BAAN65152	S-BAAN65252	S-BAAN65154

BONE LEVEL SHORT IMPLANT COUPLE ABUTMENTS



ø5.0



ø6.0



ø8.0

	2 mm	3 mm
ø5.0	S-BMCA5020	S-BMCA5030
ø6.0	S-BMCA6020	S-BMCA6030
ø8.0	S-BMCA8020	S-BMCA8030

TEMPORARY ABUTMENTS



only for ø3.3 implant

S-BMTA37

ø4.5

S-BTA45

ø5.5

S-BTA55

BONE LEVEL SHORT IMPLANT IMPRESSION COPING PICK UP



ø5.0

S-BMSICP50



ø6.0

S-BMSICP60

BONE LEVEL IMPRESSION COPING TRANSFER



only for ø3.3 implant

	HEX	NON-HEX
ø3.7	S-BMICT37L S-BMICT37S	S-BMICTN37L S-BMICTN37S
ø4.5	S-BICT45L S-BICT45S	S-BICTN45L S-BICTN45S
ø5.5	S-BICT55L S-BICT55S	S-BICTN55L S-BICTN55S
ø6.5	S-BICT65L S-BICT65S	S-BICTN65L S-BICTN65S

BONE LEVEL IMPRESSION COPING PICK UP



only for ø3.3 implant

	HEX	NON-HEX
ø3.7	S-BMICP37L S-BMICP37S	S-BMICPN37L S-BMICPN37S
ø4.5	S-BICP45L S-BICP45S	S-BICPN45L S-BICPN45S
ø5.5	S-BICP55L S-BICP55S	S-BICPN55L S-BICPN55S
ø6.5	S-BICP65L S-BICP65S	S-BICPN65L S-BICPN65S

Ti-Base Abutment



S-BMTIA05
S-BMTIA10
S-BMTIA20



S-BTIA05
S-BTIA10
S-BTIA20

BONE LEVEL SHORT IMPLANT IMPRESSION COPING TRANSFER



S-BMSICT50



S-BMSICT60

BONE LEVEL FIXTURE LAB ANALOG



S-BMFLA
only for ø3.3 implant



S-BFLA

HEALING ABUTMENTS

Bone Level Healing Abutments

				
	ø3.7 only for ø3.3 implant	ø4.5	ø5.5	ø6.5
1 mm	-	S-BHA4510	S-BHA5510	S-BHA6510
2 mm	-	S-BHA4520	S-BHA5520	S-BHA6520
2.5 mm	S-BMHA3725	-	-	-
3 mm	-	S-BHA4530	S-BHA5530	S-BHA6530
3.5 mm	S-BMHA3735	-	-	-
4 mm	-	S-BHA4540	S-BHA5540	S-BHA6540
5 mm	S-BMHA3750	S-BHA4550	S-BHA5550	S-BHA6550
7 mm	S-BMHA3770	-	-	-

Bone Level Short Implant Healing Abutments

			
	ø5.0	ø6.0	ø8.0
2 mm	S-BMHA5020	S-BMHA6020	S-BMHA8020
3 mm	S-BMHA5030	S-BMHA6030	S-BMHA8030

O-RING SYSTEM

O-RING ABUTMENTS

				
	ø3.42	ø4.5	ø2.9 only for ø3.3 implant	ø3.7 only for ø3.3 implant
0.5 mm	S-BORA05	-	S-BMORA05	-
1 mm	S-BORA10	-	-	-
2 mm	-	S-BORA20	-	S-BMORA20
3 mm	-	S-BORA30	-	-
4 mm	-	S-BORA40	-	S-BMORA40
6 mm	-	S-BORA60	-	-

RETAINER



S-ORR

O-RING LAB ANALOG



S-ORLA



GREEN RETAINER
Retention 350 Gr



TITANIUM CAP

S-ORRo

LOCATOR SYSTEM

BONE LEVEL MINI LOCATOR ABUTMENTS

	
	ø3.3 only for ø3.3 implant
1.5 mm	S-BMLA15
2.5 mm	S-BMLA25
3.5 mm	S-BMLA35
4.5 mm	S-BMLA45
5.5 mm	S-BMLA55
6.5 mm	S-BMLA65

BONE LEVEL LOCATOR ABUTMENTS

	
	ø3.86
1.5 mm	S-BLA15
2.5 mm	S-BLA25
3.5 mm	S-BLA35
4.5 mm	S-BLA45
5.5 mm	S-BLA55
6.5 mm	S-BLA65

LOCATOR MATRIX



S-LMHt



VIOLET CAP
Rigid Retention 2.7 Kg



YELLOW CAP
Extra Soft Retention 0.6 Kg



CLEAR CAP
Standard Retention 1.8 Kg



BLACK CAP
Processing



PINK CAP
Soft Retention 1.2 Kg

LOCATOR LAB ANALOG



S-BLLA Ø 4mm

Narrow Implant Accessories


NARROWIMPLANT

It is developed for use in the anterior region, where the mesio-distal distance is not sufficient, or in extremely thin crest. It can be preferred in aesthetically challenging situations.

THE NARROW IMPLANT Ø 2.9 mm IS ONLY INDICATED FOR:

Single tooth restorations:

- Central & lateral incisors in the mandible
- Lateral incisors in the maxilla



Bone Level Narrow Couple Abutment

Ø3.5

2.0mm	S-BNCA35207
3.0mm	S-BNCA35307
4.0mm	S-BNCA35407
6.0mm	S-BNCA35607



Bone Level Narrow Angled Abutment

Ø3.5

15°x2(Cuff)	S-BNAA35152
15°x3(Cuff)	S-BNAA35153



Bone Level Narrow Temporary Abutment

Ø3.5

14.85mm	S-BNTA35
---------	----------



Bone Level Narrow Healing Abutment

Ø3.5

3mm	S-BNHA3530
6mm	S-BNHA3560



Bone Level Narrow Abutment Screw

Ø1.90

6.5mm	S- BNAS
-------	---------



Bone Level Narrow Closing Screw

Ø2.60

4.5mm	S-BNCLS
-------	---------



Bone Level Narrow Impression Coping Pick Up

Ø3.5

S- BNICP37
S- BNICPGL



Bone Level Narrow Impression Coping Transfer

Ø3.5

S- BNICT37
S- BNICTGS



Bone Level Narrow Lab Analog

Ø3.5

S-BNFLA




Bone Level Narrow Ti-Base Abutment

Ø3.5



0.5mm	S-BNTiA05
1mm	S-BNTiA10
2mm	S-BNTiA20




STANDARD KIT & PROTOCOLS



SHORT KIT & PROTOCOLS





ZYGOMATIC KIT



SURGICAL DRILLS

ZYGOMA DIAMOND BURR		S-ZBURN
ZYGOMA POINT DRILL		S-PD2970
ZYGOMA TWIST DRILL 2,5 SHORT		S-TWD3565
ZYGOMA TWIST DRILL 2,5 LONG		S-TWD35100
ZYGOMA TWIST DRILL 2,9 SHORT		S-TWD2965
ZYGOMA TWIST DRILL 2,9 LONG		S-TWD29100






GUIDED KIT & PROTOCOLS



GUIDED SURGERY SYSTEM


Implantswiss Guided Surgery System is designed to provide reliable, precise and fast solutions in flapless surgery even in the complicated indications.

- Safe
- Precise
- Less Surgery Time
- Minimal Invasive



Experience

The Perfection



 02 0102 3444 446

 23 Dr. El mahrouky st. Agouza 7th floor flat 72. Giza. Egypt

 02 23 303 3934

 www.verdentegypt.com  info@verdentypt.com

 02 23 303 3926

   Verdent Egypt